OnGuard[®] 2 Closed System Transfer Device (CSTD)

Containing Chemfort[™] CSTD Reconstitution and Administration Components

P-8749-2 REV. 09/22

VIAL ADAPTORS (VA) AND CONVERTOR 13 mm

DRUG █

DILUENT/SALINE █

13 mm	20 mm	28 mm
VA 20 mm + convertor 13 mm	VA 20 mm	VA 28 mm

SYRINGE ADAPTOR (SA)

SYRINGE ADAPTOR LOCK (SAL)

LUER LOCK ADAPTOR (LLA)

BAG ADAPTOR SP (BASP)

VIAL ADAPTOR CONNECTION
Each Vial Adaptor fits a specific vial neck

A

20mm / 28mm
13mm
Remove cap

20mmVA / 28mmVA

2a

3a

CLICK!

20mm / 28mm
Connect VA spike
Once connected, do not disconnect

B

SA/SAL CONNECTION

1
Male Luer Lock (MLL) syringe

2a

Connect SA to syringe
Once connected, do not disconnect

2b

Connect SAL to syringe

C

BAG ADAPTOR SP

1
Remove cap

2

3

Insert BASP to Infusion Bag
Once connected, do not disconnect

DRUG RECONSTITUTION (powder for injection)

D

13mm / 20mm / 28mm
Remove caps

2a

CLICK!

Connect SA/SAL to VA of diluent, invert and withdraw diluent

1b

2b

3b

Connect SA/SAL to BASP; Withdraw diluent; Pinch wings to disconnect

4
5
6
7

Connect SA/SAL of diluent to VA of lyophilized vial, inject diluent; dissolve drug

DRUG TRANSFER

E

13mm / 20mm / 28mm
Remove caps

2
CLICK!

Connect SA/SAL to the VA of drug and invert to withdraw

3
Withdraw drug while vial is upwards
Discard air bubbles. Do not inject fluid upwards

4
Pinch wings on SA/SAL to disconnect

5
Re-Cap SA/SAL and VA for future use

DRUG INFUSION TO IV BAG

F

1
CLICK!

Connect SA/SAL of drug to BASP and inject drug

2
3
4
5
2-3 mL

Flush with 2-3 mL of solution

Pinch wings on SA/SAL to disconnect

ADMINISTRATION-IV DRIP

G

1a
2a

USING BASP

Open Spike Port tail and insert primed Set Spike
Once connected, do not disconnect

ADMINISTRATION- IV PUSH/BOLUS

H

1
FLL

2
LLA

Connect LLA to IV Female Luer Lock (FLL) port
Once connected, do not disconnect

3
CLICK!

5 sec IPA 70%

4
Connect SA/SAL of Drug to LLA and inject

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Explanation of symbols taken from ISO 15223-1:2021

ref: 5.1.1		Manufacturer	ref: 5.1.6		Catalog number	ref: 5.2.8		Do not use if package damaged	ref: 5.3.2		Keep away from sunlight	ref: 5.4.3		Consult instructions for use
ref: 5.1.4		Use-by date	ref: 5.2.3		Sterilized using ethylene oxide	ref: 5.2.11		Single sterile barrier system	ref: 5.3.4		Keep dry	ref: 5.6.3		Non-pyrogenic
ref: 5.1.5		Batch code	ref: 5.2.6		Do not re-sterilize	ref: 5.2.14		Single sterile barrier system with protective packaging outside	ref: 5.4.2		Do not re-use	ref: 5.7.8		Translation